

Technical sheet

MaxiBio 50

Disinfection by misting - compressed air

Agis Medical Supplies Pte Ltd

Tel: 62599788

Email: sales@disinfectant.sg



MaxiBio 50 can disinfect surfaces by air by saturating a volume with a disinfectant.

The disinfection of surfaces by air is fast, simple and very effective. The MaxiBio 50 fogger allows you to carry out a disinfection of the surfaces by air of the premises of up to 50 m³.

This mist generates in the atmosphere an ecological disinfectant with exceptional qualities :

Hydrogen Peroxide

It is ideal for the disinfection of :

- Refrigerated trucks, sanitary vehicles, ambulances
- Medical, veterinary, dental practices
- Nurseries (toys, dormitories, furniture, diapers etc ...)
- Housing after death
- Microbiological safety stations (PSM), laboratories
- Ski boots, bowling shoes
- Nursing homes, rooms of a retirement home after death or leaving a patient

Dry misting is a process that makes it possible to saturate a room with a disinfectant: all surfaces in contact with the air are treated without human or animal presence.

Technical sheet

| | |
|------------------|--------------------------------|
| Alimentation | 220-240V / 50 Hz |
| Maximum Power | 240V / 50 Hz |
| Dimensions | H 260 x l 330 x P 110 mm |
| Weight | 6 kg |
| Product flow | 10 ml / min |
| Disinfectant | Hydrogen Peroxide |
| Compressor | Piston lubrifiant non lubrifié |
| Air flow | 27 NL / min |
| Rotation speed | 2800 rpm |
| Sound level | 56 dB |
| Mode | Automatic |
| Broadcast Time | Adjustable |
| Deferred startup | 30 s to 1 min |
| Maximum Volumes | 50 m ³ |

Hydrogen Peroxide Disinfectant



Powerful disinfectant based on hydrogen peroxide and peracetic acid, non-toxic to humans and their environment.

Non- polluting solution and biodegradable.

No aldehyde, no phenol, no hypochlorite, no heavy metals.